

Art.No.	Description	Li	st price		Offer price			
	Description	SEK/NOK	DKK	EUR	SEK/NOK	DKK	EUR	
BRO-46-0001	IntelliXcap S1 Semi-Automated Screw Cap Decapper Re-capper	19995	14995	1820	Fixed	ce		
430791	Centrifuge tube, PP, 15 ml, sterile 500 pcs, bulk	2839	1987	258	605	395	53,04	
430790	Centrifuge tube, PP, 15 ml, sterile 500 pcs, rack	2866	2006	261	617	403	54,09	
430829	Centrifuge tube, PP, 50 ml, sterile 500 pcs, bulk	2951	2065	268	621	406	54,48	
430828	Centrifuge tube, PP, 50 ml, sterile 500 pcs, rack	3298	2309	300	633	414	55,53	
430921	Centrifuge tube, PP, 50 ml, sterile 500 pcs, self standing, bulk	3560	2492	324	921	602	80,78	



LIQUID HANDLING



ERGONOMIC MOTORIZED PIPET CONTROLLER

- for fast and accurate Liquid Handling

This light, ergonomic, motorized pipet controller is designed for use with glass and plastic serological pipets from 0.5 to 100 ml.

Conveniently positioned switches change operating modes and speeds quickly to handle different liquid volumes and various viscosities.

> Aspirating and dispensing speed is controlled through the finger triggers.

Incl.: 0.2 µm filter, universal power supply, li-ion battery set, and two-position charging stand.

Art.No.	Dominion		List price		Offer price			
	Description	SEK/NOK	DKK	EUR	SEK/NOK	DKK	EUR	
<u>357801</u>	Stripettor Pro Pipet Controller	4700	3290	428	2326	1520	204,00	
4487	Serological Pipets, 5 ml, paper/plastic wrap, 200 pcs	1054	737	96	290	189	25,42	
4488	Serological Pipets, 10 ml, paper/plastic wrap, 200 pcs	1072	751	97	303	198	26,59	
4489	Serological Pipets, 25 ml, paper/plastic wrap, 200 pcs	1995	1396	181	539	353	47,32	
4490	Serological Pipets, 50 ml, paper/plastic wrap, 100 pcs	2214	1550	201	571	373	50,06	

PRECISION IN EVERY DROP ♦ ♦



Starlabs RepeatOne S is perfect for routine or research applications in diagnostics, molecular biology, and environmental analysis.

- √ Up to 49 dispenses from the same sample
- \checkmark Dispensing volume 2.0 μ l 5.0 ml
- √ Easy tip insertion tip is simply inserted from below
- √ Increased chemical resistance through innovative plastic materials
- Conforms to EN ISO 8655 Piston-Operated Volumetric Apparatus

STROKE INDICATOR

BUY A DISPENSER AND GET **TIPS FOR FREE**



Art.No.	Description		List price			Offer price	e	
	Description	SEK/NOK	DKK	EUR	SEK/NOK	DKK	EUR	
G4026-1500	RepeatOne S dispenser, incl. shelf mount	4464	3125	406		enser and f tips for fre	_	
<u>S4761-0500</u>	Repeat despenser tip 5 ml sterile, 100 pcs	1735	1215	158	1475	1033	134,3	
<u>S4761-1000</u>	Repeat despenser tip 10 ml sterile, 100 pcs	1735	1215	158	1475	1033	134,3	
<u>S4761-2500</u>	Repeat despenser tip 25 ml sterile, 100 pcs	2506	1754	228	2130	1491	193,8	
<u>S4761-5000</u>	Repeat despenser tip 50 ml sterile, 100 pcs	2575	1803	234	2189	1533	198,9	



CONSUMABLES

CHOOSE YOUR STARLAB STARGUARD GLOVES

Best suited to the job in hand

COMFORT - COMFORTABLE AND STRONG

- √ PPE Cat. III Complex Design
- √ 245 mm long
- √ Textured and snug fitting fingers
- √ Made from soft nitrile

Comfort - Art.No. SG-C-XS, -S, -M, -L, -XL

PROTECT - OUTSTANDING BREAK-THROUGH **TIMES**

- √ n-Heptane > 60 mins
- √ Sodium hydroxide (40 %) > 480 min
- Formaldehyde (37 %) > 240 min
- Hydrogen peroxide (30 %) > 30 min

Protect - Art.No. SG-P-XS, -S, -M, -L, -XL



FROM 04 SEK/NOK 0.73 DKK

ALL STARGUARD GLOVES ARE MANUFACTURED TO THE **HIGHEST** STANDARDS AND ARE TESTED AND CERTIFIED BY A THIRD-PARTY NOTIFIED BODY.

Ask for more information, samples and other variants.



DID YOU KNOW THAT.

Two of Starlabs gloves have an Acceptable Quality Level (AQL) of 0.65.

For EN ISO 374-5 for the protection against bacteria and fungi, all gloves also protect against viruses, and they are all in compliance with EN safety standards.

For EN ISO 374-1, all gloves are classified as Type B for the requirements for chemical risks.

All gloves have been tested against a variety of chemicals to EN 16523-1:2015+A1:2018 standards.

FLUIDX STORAGE TUBES

Superior tracking and automation for sample audits

FluidX sample tubes make it easy to track and share data between users, labs, and locations, with strong for automation.

The ranges of FluidX tubes includes:

- √ External and internal thread tubes
- 2D-coded tubes
- √ Dual-coded tubes
- $\sqrt{}$ Format 24, 48 and 96 wells
- √ Working volume 0.26 -7.6 ml



Art.No.	Description		List price			Offer price		
	Description	SEK/NOK	DKK	EUR	SEK/NOK	DKK	EUR	
BRO-65-7642	FluidX 48 tubes, 1.9 ml, Ext. Thread, Tri-coded, 2D, 10 racks	9323	6526	848	8599	6019	781	
BRO-68-0703-11	FluidX 96 Tubes, 0.5 ml, Ext. Thread, Tri-coded, capped, 10 racks	16449	11500	1495	14523	10166	1320	
BRO-68-1003-11	FluidX 96 Tubes, 1.0 ml, Ext. Thread, Tri-coded, capped, 10 racks	15740	11019	1431	14523	10166	1320	
BRO-68-0701-11	FluidX 96 tubes, 0.5 ml, Ext.Thread, Dual-coded, capped, 10 racks	11074	7752	1007	10214	7150	929	



INSTRUMENTS



• Sample capacity: 10 - 100 µl

Temp. Range: 4°C~105°C min. resolution: 0.1°C

• Hot-lid temp. range: 30°C - 110°C

• Mobile App Compatibility: iOS / Android

• Temp. Accuracy: ≤ 0.2°C

• Memory: 2.000, USB flash support

• **CE-IVD** certified version is also available



GENE EXPLORER THERMAL CYCLER Intuitive software, flexible ramp rate

- √ Intelligent Operation A new generation of CPU, 8-inch responsive touch screen, and userfriendly UI design makes your operation simple and faster.
- √ Temperature Accuracy Newly-designed block and tailor-made peltier to ensure exceptional thermal performance and optimize your PCR temperature with good accuracy.
- √ Mobile App Sync Innovative mobile app and WiFi capability helps to monitor your PCR status with remote control.

SEMI-AUTOMATED SHEET HEAT SEALER Compatible with a wide range of seals and plates

Plate and seal are placed on the holder, the "Operate" button pressed and the drawer automatically closes. The sealing process is controlled by an electric mechanism and this unique feature eliminates operator variation, giving a consistent seal every time and

reducing the damage caused to microplates by over-sealing.

- 1. Set time and temperature
- 2. Place a microplate in the open drawer
- 3. Place a heat seal on top of the plate and press the button





20 % DISCOUNT



Standard plate carrier included.



CORIOLIS COMPACT

Portable dry air sampler

The Coriolis Compact is a dry cyclonic collector intended for microbial air monitoring. Light and compact, it can easily be transported and positioned on strategic places for indoor and outdoor sampling.

- 1. Collection time: From 1 min 8 hours
- 2. Flow: 50 L/min
- 3. Size of particles collected: 500 nm 10 µm
- 4. Temp. range: +5°C +45°C

SEK/NOK 68.500 DKK 8.900 €

Art.No. P002055-CORCO-A.0



THE LOCAL ALTERNATIVE IN THE NORDICS FOR EVERYONE WHO WORKS IN LABORATORIES – BIG ENOUGH TO BE ABLE, SMALL ENOUGH TO CARE

INSTRUMENTS

MINI NUTATING MIXER

Combine Shaking and Rocking

Combines the motions of an orbital shaker and a rocker to produce three dimensional action. The benefit of this type of motion is that it is both gentle and thorough.

√ **Motion:** 3-dimensional rocking rotating

Speed: 20 rpm

√ Temp. range: 5°C - 40°C									
Art.No.	Description	List price Offer		fer pric	er price				
	Description	SEK/NOK	DKK	EUR	SEK/NOK	DKK	EUR		
<u>S0600-230V</u>	Labnet MiniMix Nutating 3-D shaker with dimpled rubber mat	5845	4090	532	5261	3681	478,80		
E2110-BBL	Labnet Blot box, 11.7 x 8.9 cm, 6 -10 ml capacity for larger mini protein gels	1605	1124	146	1445	1012	131,40		
E2110-BBS	Labnet Blot box, 9.1 x 6.6 cm, 3-50 ml capacity for mini protein gels	1458	1021	133	1313	919	119,70		



SHAKEN - NOT STIRRED

Similar shakers with different movements

See-saw rocking or 3D mixing, both shakers have a fixed tilt angle and provide two shaking modes: continuous or timer. Both shakers are supplied with a platform with bungee cords.

GET SHAKING

- √ Speed range 10-70 rpm
- √ Movement 3D Mixing / See Saw Rocking
- √ LCD display for speed and time
- √ 2 modes continuos and timer
- √ Alternative platforms available



Art.No.	Domintion	L	ist price		C	Offer price	
	Description	SEK/NOK	DKK	EUR	SEK/NOK	DKK	EUR
N2400-8010	Rocking Shaker	18900	13230	1718	14595	10200	1327
N2400-8020	3D Rocking Shaker	18900	13230	1718	14595	10200	1327



DID YOU KNOW THAT....

Whether you are performing routine analyses or specialized procedures, we offer a wide range of essential laboratory instruments designed to support all your daily lab needs. Learn more: swab.se/en/Products/Instruments







ANYTIME ANYWHERE - EASY SAMPLE MEASUREMENTS WITH A PORTABLE **pH METER**

Endress + Hauser 🖽

Whether you are in the lab or out measuring in streams, waterways, fish farms or wastewater treatment plants, the handheld mobile CML18 makes it easy for you.

Set incl.

- √ CML18 mobile device measurement
- √ CPL51E Memosens 2.0 pH sensor
- √ CML18 protective cover for the mobile device
- √ 1,5m CYK20 flexible data cable
- √ Fieldcase
- √ pH buffer 4 and 7 for all industries


CHOOSE THE 3D CELL CULTURE MODEL 3 That works best for your research

Whether you are starting 3D cell culture, scaling up, or shifting to high-throughput screening, Corning products help you create more in vivo-like, predictive models quickly and efficiently.

→ **Spheroids** can be scaled up using specialized multiwell growth plates. As cell clusters, they form layers. External cells have the most contact with media, drugs, and nutrients.

→ **Organoids** grow as mini-organs, maintaining cell functionality or characteristics similar to in vivo. They are more complex to grow but behave more naturally than 2D cultures.

→ **Tissue Models** have been developed for a variety of issues to display micro anatomy, mimic complex tissues and organ function, and shed light on cell-to-cell interactions.









WE SEE A GROWING DEMAND FOR GREENER SOLUTIONS FROM OUR CUSTOMERS

To meet this need we are proud to collaborate with a supplier who has a prioritized sustainability for many years, with long-term corporate planning and a strong focus on producing durable, highquality products.

Starlab has already succeeded in saving plastic with their TipOne pipette tip system, and the journey to make pipetting more sustainable still continues.



REDUCE AND REUSE

- 🏲 Material TipOne racks, refills and stack racks are made entirely from recyclable polypropylene and contain no paper labels, so recycling is quick and easy.
- Green energy TipOne is manufactured in a factory powered entirely by wind.
- Recycled plastic All clear plastic packaging on the refills (bases, lids and spacers) are manufactured from 100% recycled platic.
- Plastic saved Up to 63%, when using a sterile refill compared to using the racked equivalent.
- Carboard boxes Pack and case cartons for TipOne are made from 100 % recycled cardboard.



SCAN THE QR **CODE FOR MORE INFORMATION** ABOUT STARLABS **GREEN MISSION**





Reference: https://www.starlabgroup.com/en/sustainability.html



BIOPROCESS



<u>Cleargreen® Pure</u>

Tubing for Advanced Therapies

Setting a New Standard in PVC Tubing for Blood Collection and Cell and Gene Therapy Processes.

Cleargreen Pure is a sustainable, plant derived vinyl material that is DEHP free, gamma stable and compatible with any aseptic welder or sealer. It has been extensively studied for its superior toxicological and physiochemical profile.

WHY CHOOSE CLEARGREEN PURE

- √ DEHP free
- √ Animal Derived Ingredients Free (ADIF)
- √ Gamma stability
- √ Rigorously studied
- √ Rigorous testing
- √ Sustainability

Contact us for more information about Sanisure's Cleargreen Pure Tubing for Advanced Therapies.







CONTACT INFORMATION

SWEDEN

PRICES
ON REQUEST

Phone: +46 40 51 00 00 Web: www.swab.se Email: info@swab.se

NORWAY

Phone: +47 22 22 87 87 Web: www.swab.no Email: info@swab.no

Opening hours

Monday - Friday: 8.00-17.00

DENMARK

Phone: +45 86 20 16 16 Web: www.swab.dk Email: info@swab.dk

FINLAND

Phone: +358 9 315 45 620 **Web:** www.swab.fi **Email:** info@swab.fi





*All prices are net and valid until December 25 and cannot be combined with other agreements or discounts. Quote discount code HK25 when ordering. General terms and conditions of sale apply, see our <u>website</u>. We reserve the right for printing errors.